APTT-EA liquid

Diagnostic reagent for quantitative partial Thromboplastin Time (APTT) determination in human plasma on coagulation instruments (ellagic acid activated)

REF	Cont.	Product	Description
C01130	5x 8 mL	APTT-EA liquid	
C01130SV	1x 8 mL	APTT-EA liquid	yellowish turbid liquid
C14130	2x 5 mL + 2x 5 mL	APTT-EA liquid 0.020 M CaCl ₂	clear colourless liquid

Additionally offered:

C03101	5x 10 mL	0.020 M Calcium Chloride
T955100	10x 1 mL	Normal Coagulation Control Plasma
T955101	10x 1 mL	Abnormal Coagulation Control Plasma

For professional in vitro diagnostic use only.

GENERAL INFORMATION

Method:	Clotting Assay
Temperature:	37°C
Sample:	Sodium Citrate Plasma

Number of tests:		
C01130	5x 80 Tests 5x 160 Tests 5x 320 Tests	Manual Method Micro Method (DiaCheck CP4) Quarter Method (DiaCheck C1, C2, C4)
C01130SV	80 Tests 160 Tests 320 Tests	Manual Method Micro Method (DiaCheck CP4) Quarter Method (DiaCheck C1, C2, C4)
C14130	100 Tests 200 Tests 400 Tests	Manual Method Micro Method (DiaCheck CP4) Quarter Method (DiaCheck C1, C2, C4)

REAGENT COMPOSITION

- Components:
- Phospholipid cephaloplastin reagent derived from rabbit brain
 - Ellagic acid

REAGENT PREPARATION

The reagent is ready to use

MATERIALS REQUIRED BUT NOT PROVIDED

- Calcium Chloride (0.020 M, REF: C03101)
- A mechanical or photo-optical Coagulation system, either manual, semi-automated or fully automated.
- Suitable timer.
- Parafilm or equivalent.
- Precision pipette, 0.1 mL
- Control Plasmas.

REAGENT STABILITY AND STORAGE

Conditions:	close immediately after use
Storage:	at 2-8°C
Stability:	up to the expiration date, once opened 3 months at 2-8° C, 1 week at 18-25° C and 2 days at 37° C

Avoid prolonged heating! Do not freeze!

SAMPLE COLLECTION AND PREPARATION

- Obtain venous blood by clean venipuncture.
- Immediately mix 9 parts of freshly collected blood with 1 part tri-sodium citrate (0.11 mol/L, 3.2%).
- Centrifuge immediately for 15 minutes at approximately 1500 g (3000 rpm).
- Transfer the plasma into a clean test tube. It should be ensured that the plasma is free from platelets (PPP)
- Plasma has to be tested preferably immediately. If the specimen are held at 22 – 24°C then they may be tested within 2 hours and if the specimen is held at 2 – 4°C then it may be tested within 3 hours.

MANUAL TEST PROCEDURE

Bring reagents and sample to 37°C

Aspirate from reagent vial enough reagent for immediate testing requirements in a thoroughly clean and dry test tube. (Plastic tubes are preferred) Pre-warm reagent at 37°C

Pipette into test tubes	Sample	Reference Plasma
Sample (prewarmed at 37°C)	100 µl	-
Reference Plasma	-	100 µl
APTT Reagent	100 µl	100 µl

Shake tube briefly to mix the reagent and plasma; place the tube at 37°C for 3 to 5 minutes. Then forcibly add:

Calcium Chloride (pre-warmed at 37°C)	100 µl	100 µl
---------------------------------------	--------	--------

Simultaneously start a stopwatch and incubate for 20 seconds at 37°C. Following incubation, gently tilt back and forth until a gel clot forms. Stop the watch and record time in seconds.

For micro method: Use ½ of each volume
For quarter method: Use ¼ of each volume

Great care must be taken to minimize variations which may occur by seemingly insignificant factors:

Calibration Curve Method (For determination of heparin concentration)

- Dilute heparin (as used for treatment) with physiological saline to a concentration of 10 U/mL.
- Mix 0.2 ml of 10 U/ mL diluted heparin with 1.8 mL of Reference Plasma to give a heparin standard of 1 U/mL concentration.
- Dilute the heparin standard as prepared above (1 U/mL) with Reference Plasma as follows:

Test tube no.	Heparin Std (1 U/ mL) in mL	Reference Plasma in mL	Heparin conc. (U/ mL)
1	0.5	-	1
2	0.4	0.1	0.8
3	0.3	0.2	0.6
4	0.2	0.3	0.4
5	0.1	0.4	0.2
6	0.1	0.9	0.1
7	-	0.5	0.0

CALCULATION

$$R = \frac{\text{APTT of patient plasma (in seconds)}}{\text{APTT of Reference Plasma (in seconds)}}$$

+ Alternatively a patient plasma pool from at least 20 normal healthy individuals can be used.

Calibration Curve:

Heparin Concentration in the test sample can be directly obtained from the calibration curve by interpolating the test plasma clotting time against the heparin concentration in U/mL.

REFERENCE RANGE*

22 – 34 seconds (at 3 minutes activation time)

* It is recommended that each laboratory establish its own normal range.

TEST PRINCIPLE

Cephaloplastin activates the coagulation factors of the intrinsic pathway of the coagulation mechanism in the presence of calcium ions. APTT is prolonged by a deficiency of one or more of the clotting factors of the intrinsic pathway and in the presence of coagulation inhibitors like heparin

QUALITY CONTROL

All control plasmas with APTT values determined by this method can be used.

We recommend:

REF	Cont.	
T955100	10x 1 mL	Normal Coagulation Control Plasma
T955101	10x 1 mL	Abnormal Coagulation Control Plasma

CALIBRATION

For the assay with Calibration Curve Method Heparin as used for Treatment is required.

All reference plasmas with APTT values determined by this method can be used. We recommend:

REF	Cont.	
T955105	10x 1 mL	Normal Coag. Reference Plasma

WARNINGS AND PRECAUTIONS

- Due to inter and intra laboratory variations users must establish their own normal population range as well as normal and abnormal range.
- It is recommended that controls with known factor activity should be run simultaneously with each test series routinely.
- Incorrect mixture of blood and tri-sodium citrate, insufficient pre-warming of plasma and reagent, contaminated reagents, glassware etc. are potential source of errors.
- Incorrect dilutions of heparin are also a potential source of errors.
- Oxalated plasma may induce prolonged clotting times
- Clotting time of patients on anticoagulant therapy depend upon the type and dosage of anticoagulant and also the time lag between the specimen collected and the last dose.
- Abnormalities of coagulation factor VII, factor XIII and platelets are not detected by this method.
- For automated equipment it is strongly recommended that the equipment manufactures methodology is strictly adhered to.
- In heparin monitoring time of collection of blood sample is important since the in-vivo half-life of heparin is approximately 1.5 hours. When it is administered intravenously it has an immediate anti-coagulant effect but its efficacy decreases rapidly with time.
- Platelet factor IV, a heparin-neutralizing factor can be released due to platelet aggregation or damage. In order to prevent this phenomenon in vitro the specimen should be collected with a minimum of trauma.
- Decrease in APTT time is observed in males under estrogen therapy and oral contraceptive administration in females.
- Glassware and cuvettes used in the test must be scrupulously clean and free from even traces of acids/alkalis or detergents.

WASTE MANAGEMENT

Please refer to local legal requirements.

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