

URICHA - Urine Chemistry Analyser

Features

- Easy to use, handheld low power analyzer.
- Finger touch operation makes easy to interface.
- Large memory to save 250 Results with date.
- Convenient computer interface via USB to back up the results.
- Single input voltage 12 VDC.
- Large display for clear view.
- Results for all 10 parameters are conveniently displayed on one page.
- Alphanumeric Input possible
- Instrument can be connected to Laboratory information systems and also to PC using USB/RS232 interfaces & software.
- Built in thermal printer with 24 characters.