



DiaMetra

care for quality

Anti-CP 4 IgG/IgA

Determination of IgG /IgA antibodies
against Linear Citrullinated **Peptides 4**

Based on the latest generation of linear citrullinated peptides

- **Improved sensitivity : 80,0%**
- **Specificity : 97,0%**



Anti CP 4 IgG/IgA

Determination of IgG /IgA antibodies against Linear Citrullinated Peptides 4

The antigens used in anti CP4 kit is an optimized combination of four new different viral linear peptides. The kit detects along with anti IgG, also antibodies anti IgA in RA patients with absence of IgG antibodies, resulting in an better sensitivity than the competing products.

Correlation anti CP4 IgG vs anti CP4 IgA

Samples NHS # 30
 RA # 50

Results

NHS	IgA +	IgA -	fisher test
IgG +	0	1	
IgG -	1	28	

RA	IgA +	IgA -	P<0,005
IgG +	21	15	
IgG -	4	10	

Code	Description
* DKO 149	Anti CP IgG
* DKO 182	Anti-CP 4 IgG/IgA

*Patented and licensed by TOSCANA BIOMARKERS s.r.l. Via Fiorentina n1, Siena/It.
PATENT : WO2011/061720 - EP09176776.4 - PCT IB2010/055347

Diametra Srl - Head Office

20090 Segrate - Milano - Italy
Via Garibaldi, 18
Ph +39.02.21.39.184 - Fax +39.02.21.33.354
V.A.T. - IT 12585340156

www.diametra.com
info@diametra.com